

Measuring the chemical potential of a two-dimensional hole gas across the metal-insulator transition

An application using Oxford Instruments Superconductivity's HelioxTL, top loading single shot ^3He refrigerator

To date, the vast majority of experimental findings have been obtained through transport measurements. However, it was expected that the existence of a true zero temperature phase transition would reveal itself in other properties of the system such as the magnetisation and electronic compressibility, $\partial n/\partial\mu$, where μ is the chemical potential and n is the carrier density.

For example, a transition to a Mott insulator would result in a vanishing compressibility at the transition, in contrast with a crossover to an Anderson insulator where the compressibility is expected to remain nearly unaffected. It was therefore clear that determining the chemical potential as a function of the carrier density and hence the compressibility across the transition would help to unravel the microscopic origin of this phase transition.

Recent experiments carried out by Drs Ilani, Yacoby, Mahalu and Shtrikman at the Weizmann Institute of Science, Israel, have measured the dependence of the compressibility of a 2D hole gas (2DHG) as a function of its density as it crosses through the metal-insulator transition. The experiments required a very low base temperature, 220 mK, to be held for up to five days in a relatively large sample space. This was achieved successfully using Oxford Instruments Superconductivity's HelioxTL 38 mm ^3He insert - although even this was being stretched to its limit, the hold time being four times its design specification. The measurements were performed using several

single electron transistors (SETs) situated directly above the 2DHG. The transistors act as extremely sensitive voltage probes that measure locally the chemical potential of the 2DHG. Monitoring different transistors across




Fig. 1: Density dependence of the resistivity: Traces from top to bottom on the right correspond to $T=600, 550, 500, 450, 400, 300$ and 250 mK. The low-T crossing point is at $n_c = 2.1 \text{ cm}^{-2}$




Fig. 2: A SEM micrograph of the SET on a schematic device structure, together with the measurement circuit.

The possibility that the zero temperature ground state of a 2D system might be metallic has been raised by experiments carried out six years ago by S. Kravchenko and co-workers in silicon two dimensional electron systems. These have revealed a metal-insulator transition in two-dimensions (2D) and clearly show a metallic like behaviour above a certain critical density. Since then, this phenomenon has been observed in many different material systems such as holes in GaAs and in SiGe, challenging the scaling theory of localisation, which predicts only an insulating phase in 2D.

Contact:

Dr Amir Yacoby, Weizmann Institute, Israel
E-mail: yacoby@wis.weizmann.ac.il



the same sample also allows determination of the spatial variations in the compressibility. Simultaneous transport measurements were conducted to ensure a precise determination of the critical density within the same sample.

In the metallic phase the measurements clearly show that the compressibility is negative. This result agrees quantitatively with previous macroscopic measurements and the predictions of the mean field Hartree Fock theory. However, at the transition into the insulating phase the compressibility becomes positive and its magnitude decreases with increasing penetration into the insulating phase. Furthermore, measurements done with different transistors on the same sample show that, although the behaviour in the metallic phase is homogeneous in space, the behaviour in the insulating phase differs from point to point across the sample, suggesting that the insulating phase is spatially inhomogeneous.

These measurements, showing a clear difference between metallic- and insulator-like behaviour, therefore support the existence of a metal-insulator transition in the two-dimensional system.




Fig. 3: The absolute value of the inverse compressibility, $\partial\mu/\partial n$, collected from five SETs on three different Hall bars from two different wafers. In the metallic regime the compressibility is negative and all the SETs show the same dependence, suggesting a homogeneous phase. In the insulating regime, the compressibility is positive and vanishes. Hence the inverse compressibility diverges.

Inset: Results from two SETs on the same device demonstrating the spatial dependence of $\partial\mu/\partial n$ on the insulating side.

Oxford Instruments Superconductivity

UK
Tubney Woods, Abingdon
Oxfordshire OX13 5QX
Tel: +44 (0)1865 393 311
Fax: +44 (0)1865 393 333

Germany
Otto-von Guericke Ring 10
D-65205 Wiesbaden
Tel: +49 6122 937 171
Fax: +49 6122 937 175

Italy
Via Leone Tolstoj 86
20098 San Giuliano
Milanese
Milan
Tel: +39 02 98 2531
Fax: +39 02 98 24407

Japan
Haseman Building
2-11-6 Tomioka, Koto-ku
Tokyo 135-0047
Tel: +81 03 5245 3261
Fax: +81 03 5245 4472

Spain
Avda Mata Piñonera, 2
28700 San Sebastian de los Reyes
Madrid
Tel: +34 91 659 0740
Fax: +34 91 654 6794

U.S.A.
130A Baker Ave. Ext.
Concord, MA 01742-2121
Tel: +1 978 369 9933
Fax: +1 978 369 6616

E-mail:
superconductivity@oxinst.co.uk

Visit our web site at
www.oxford-instruments.com

This publication is the copyright of Oxford Instruments Superconductivity Limited and provides outline information only which (unless agreed by the company in writing) may not be used, applied or reproduced for any purpose or form part of any order or contract or be regarded as a representation relating to the products or services concerned. Oxford Instruments' policy is one of continued improvement. The company reserves the right to alter, without notice, the specification, design or conditions of supply of any product or service. Oxford Instruments acknowledges all trade marks and registrations

© Oxford Instruments Superconductivity Limited, 2002. All rights reserved.

Reference No: OI 600/0502

UK sites of Oxford Instruments Superconductivity Limited, operate a Quality Management System approved to the requirements of BS EN ISO 9001.

