



The detection of nickel in synthetic diamond

An application using Oxford Instruments ED2000 EDXRF spectrophotometer and MicrostatN cryostat

The Certificates Department of the Diamond High Council in Antwerp, Belgium, only issues certificates for polished, natural gem diamonds. The detection of synthetic diamonds is hence crucial.

In nature, diamond growth requires high pressure and high temperature. In the fifties scientists were, for the first time, able to imitate these growth conditions using large presses. These days an alloy of iron-nickel, which serves as a solvent-catalyst for diamond growth, is often used. The nickel is included in the diamond as atomic as well as micrometer-sized impurity. Both types of impurities can be detected by the Oxford Instruments ED2000 EDXRF spectrophotometer (see figure 1). As a substitutional atomic impurity, nickel has an easily detectable line in the photoluminescence spectrum. Strong quenching of the luminescence makes the spectrum not easy to detect at room temperature, but the signal is very strong at 77 K.

When a section of a synthetic diamond grown in a Bars press is made, yellow, colourless and brown zones (see figure 2) can be observed. This colour variation corresponds to different impurity concentrations in different crystallographic growth sectors. The yellow colour is related to single substitutional nitrogen defects in the diamond lattice, while it is suspected that the nickel causes the brown colour¹.

Luminescence is excited by a 514.5 nm Argon-ion laser and detected by a Renishaw RM2000 system with Oxford Instruments **MicrostatN** cryostat. The large optical window and long thermal stability are excellent for line scans on macroscopic samples. In the section of synthetic diamond, the intensity of the Zero Phonon Line at 493.9 nm of the nickel related luminescence is monitored in three different line scans and it appears that its intensity correlates with the brown colour (figure 3). The evidence is indirect, but the measurement supports the view that nickel causes the brown colour in the synthetic diamonds, grown from Fe-Ni catalyst-solvent.




Figure 1: The EDXRF spectrum of the synthetic diamond. The spectrum indicates the presence of nickel (Ni) and iron (Fe) because the diamond was grown in a catalyst-solvent containing iron and nickel. The Argon line (Ar) is due to the argon in the air. The presence of copper (Cu) and zinc (Zn) is due to minute traces from brass clamps in which the diamond was set prior to the EDXRF measurement.

The Certificates Department of the Diamond High Council in Antwerp, Belgium, issues certificates for natural gem stones only. The Department distinguishes between natural and synthetic stones using spectrophotometry at 77 K, which detects nickel impurities from catalysts used in the production of synthetic gems.

Reference

¹ Lawson S.C., Kanda H., Sekita M., *Phil. Mag. B*, (68), 39-46 (1993)

Contact:

Mr. Filip De Weerd, Diamond High Council (HRD), Antwerp, Belgium






Figure 2: A section of a synthetic diamond. The zones with different colours are located in different growth sectors.





Figure 3: Three line scans with a stepsize of 5 μm , and between each line scan there is a space of 500 μm . The intensity of the nickel related luminescence line seems to correlate well with the brown colour in the diamond.

Oxford Instruments Superconductivity

UK
Tubney Woods, Abingdon
Oxfordshire OX13 5QX
Tel: +44 (0)1865 393 200
Fax: +44 (0)1865 393 333

China
Rm. 14-F, No.1 Plaza
No.800 Nanjing East Road
Shanghai 200001
Tel: +86 21 63608530/1/2/3
Fax: +86 21 63608535

Germany
Otto-von Guericke Ring 10
D-65205 Wiesbaden
Tel: +49 6122 937 171
Fax: +49 6122 937 175

Italy
Via Leone Tolstoj 86
20098 San Giuliano
Milanese
Milan
Tel: +39 02 98 2215
Fax: +39 02 98 2414

Japan
Hasegan Building
2-11-6 Tomioka, Koto-ku
Tokyo 135-0047
Tel: +81 03 5245 3261
Fax: +81 03 5245 4472

Spain
Avda Mata Piñonera, 2
28700 San Sebastian de los Reyes
Madrid
Tel: +34 91 659 0740
Fax: +34 91 654 6794

U.S.A.
130A Baker Ave. Ext.
Concord, MA 01742-2121
Tel: +1 978 369 9933
Fax: +1 978 369 6616

E-mail:
superconductivity@oxinst.co.uk

Visit our web site at
www.oxford-instruments.com



This publication is the copyright of Oxford Instruments Superconductivity Limited and provides outline information only which (unless agreed by the company in writing) may not be used, applied or reproduced for any purpose or form part of any order or contract or be regarded as a representation relating to the products or services concerned. Oxford Instruments' policy is one of continued improvement. The company reserves the right to alter, without notice, the specification, design or conditions of supply of any product or service. Oxford Instruments acknowledges all trade marks and registrations

© Oxford Instruments Superconductivity Limited, 2004. All rights reserved.

Reference No: DF OI650 0204

UK sites of Oxford Instruments Superconductivity Limited, operate a Quality Management System approved to the requirements of BS EN ISO 9001.



Certificate No Q4118

