

UHV Ultrastat mounted on stand

## Features

- Use of UHV materials and fabrication techniques for total UHV compatibility
- Variable temperature operation between 3.8 K and 500 K
- Small, flexible, highly efficient cooling platform
- Simple operation with fully automated control of helium
- Operates in any direction
- Low helium consumption

## Benefits

- UHV Clean environment
- Simple to install and operate
- Low cost of ownership

## UHV Ultrastat

### Ultra High Vacuum Cryogenic Environment

#### Introduction

The UHV Ultrastat is a continuous helium flow cryostat for electrical and surface property sample measurements in a ultra high vacuum environment.


#### Description

The UHV Ultrastat is constructed from stainless steel and OFHC copper. Joints are manufactured using minimum amounts of zinc/cadmium free solder. Mounting is via a DN40CF (70 mm outer diameter) conflat flange.

A heater and rhodium-iron sensor are fitted to the cold tip. The system can be baked to 500 K. The sample may be mounted directly onto the heat exchanger for optimum thermal contact, alternatively the cryostat may be used as a 'cold finger' with a copper braid used for connection between the sample and the face of the cold tip.

#### Applications

- Measurement of electrical and surface properties of sample in UHV environments.
- Materials characterisation of nano-devices, suitable for testing and examination of semiconductors.



UHV Ultrastat can be supplied with alternate DN-CF mounting flanges and alternate length inserts.



## Standard Specification

UHV Ultrastat		
<b>UHVULTRASTAT</b>	Temperature Range	3.8 K to 500 K
	Temperature Stability	+/- 0.1 K
	Maximum Bakeout Temperature	500 K
	Cooldown Time (from ambient)	40 mins
	Cooldown Helium Consumption	1 litre
	Helium Consumption	1.1 litres hr <sup>-1</sup> at 4.2 K 0.1 litres hr <sup>-1</sup> at 10 K <0.5 litres hr <sup>-1</sup> at 150 K
	Maximum Cooling Power	1 W at 4.2 K
	Weight	1 kg

## System Components

Part Number	Name	Description
<b>UHVULTRA</b>	UHV Continuous Flow Cryostat	
<b>LLT700/13</b>	Low loss technology flexible transfer tube	Designed for use with continuous flow cryostats. Lightweight and flexible design with 180° configuration for vertical cryostat entry. 1.1 metre flexible length section. Gas cooled radiation shield, manually operated needle valve, upgradable to automated type.
<b>GF3</b>	Gas flow pump	
<b>VC31</b>	Gas flow controller	Includes vacuum gauge, helium flow meter and needle valve.
<b>ITC502</b>	Digital temperature controller	Single channel IEEE488 and RS232 interfaces. Standard calibration curves built in, auto-ranging, auto-PID setting, ability to control needle valves.
<b>SV12</b>	Storage vessel	Top-fitting for LLT and TTL type transfer tubes which have standard 12 mm dewar leg, including control valve and bladder fitting.
<b>SK</b>	Spares kit	O-rings, indium wire, screws etc.

## Optional Items

Part Number	Name	Description
<b>LLT700AUTO</b>	Upgrade from LLT700/13	Allows automatic operation of transfer tube needle valve, requires ITC controller, includes cable and monitor assembly.
<b>HD30</b>	Liquid helium dewar	Top fitting with trolley.

This publication is the copyright of Oxford Instruments Superconductivity Limited and provides outline information only which (unless agreed by the company in writing) may not be used, applied or reproduced for any purpose or form part of any order or contract or be regarded as a representation relating to the products or services concerned. Oxford Instruments' policy is one of continued improvement. The company reserves the right to alter, without notice, the specification, design or conditions of supply of any product or service. Oxford Instruments acknowledges all trade marks and registrations

© Oxford Instruments Superconductivity Limited, 2003. All rights reserved.

Reference No: OI622 02/03

UK sites of Oxford Instruments Superconductivity Limited, operate a Quality Management System approved to the requirements of BS EN ISO 9001.



Certificate No Q4118

### Oxford Instruments Superconductivity

UK  
Tubney Woods, Abingdon  
Oxfordshire OX13 5QX  
Tel: +44 (0)1865 393 200  
Fax: +44 (0)1865 393 333

China  
Room 8418a, Building A  
Jiahua Business Center  
808 Hongqiao Road  
Shanghai 200030  
Tel: +86 21 64867993  
Fax: +86 21 64861781

Germany  
Otto-von Guericke Ring 10  
D-65205 Wiesbaden  
Tel: +49 6122 937 171  
Fax: +49 6122 937 175

Italy  
Via Leone Tolstoj 86  
20098 San Giuliano  
Milanese  
Milan  
Tel: +39 02 98 2531  
Fax: +39 02 98 24407

Japan  
Haseman Building  
2-11-6 Tomioka, Koto-ku  
Tokyo 135-0047  
Tel: +81 03 5245 3261  
Fax: +81 03 5245 4472

Spain  
Avda Mata Piñonera, 2  
28700 San Sebastian de los Reyes  
Madrid  
Tel: +34 91 659 0740  
Fax: +34 91 654 6794

U.S.A.  
130A Baker Ave. Ext.  
Concord, MA 01742-2121  
Tel: +1 978 369 9933  
Fax: +1 978 369 6616

E-mail:  
superconductivity@oxinst.co.uk

Visit our web site at  
[www.oxford-instruments.com](http://www.oxford-instruments.com)

