

**X-Max<sup>®</sup>**


The largest area SDD

# SIZE IT MATTERS

20mm<sup>2</sup>

50mm<sup>2</sup>

80mm<sup>2</sup>



For accuracy / For nano-analysis  
For throughput / For productivity

**OXFORD**  
INSTRUMENTS

*The Business of Science<sup>®</sup>*




# X-Max

## Large area SDD

### Analysis with a large area analytical SDD benefits from:

- Simultaneous imaging and analysis without compromise
- Large active area with superb analytical resolution
- Practical nanoanalysis at productive count rates
  - High count rates, even at low kV
  - High count rates, even at small spot sizes
- The correct results

A traditional 10mm<sup>2</sup> detector has a small capture angle = fewer counts  
A modern large area detector has a large capture angle = many counts




20mm<sup>2</sup>

50mm<sup>2</sup>

80mm<sup>2</sup>

### X-Max delivers:

- Up to 1,000% greater productivity
- Real 200,000 cps analysis
- Usable nanoscale analysis
- The same excellent resolution and results at all sensor sizes



— Traditional SDD  
— X-Max SDD

**X-Max** The largest area SDD

# X-Max

Productivity and accuracy combined


20mm<sup>2</sup>, 50mm<sup>2</sup>, and 80mm<sup>2</sup> SDDs with identical world-beating analytical performance

Now you can have count rate, imaging, and analytical performance all at the same time:

- Up to 80mm<sup>2</sup> active area
- Count rates > 500,000 cps
- Throughput > 200,000 cps

Maximum counts from larger detectors with increased solid angle:

- Count for less time - increase productivity
- Or count for the same time and get more precise results
- Or gather data at much lower beam currents meaning:
  - No sample damage
  - Longer filament tip lifetimes
  - Reduced contamination
  - Better spatial resolution
- Practical nanoanalysis with productive count rates at low accelerating voltages



**X-Max** analytical SDDs offer the best solution for productive analysis

## X-Max features


- Unique large area SDD sensors
- Discrete PentaFET® design - accuracy at all sizes
- Sensor enclosed in self-contained vacuum assembly - no oxygen X-ray absorption
- Unique electron trap for maximised solid angle
- Only one pulse processing channel required
- Tube diameter no bigger than a 10mm<sup>2</sup>
- Motorised slide as standard

**SIZE**  
IT MATTERS  
IT MATTERS

# X-Max

## Large area SDD

### Why size matters



Discover the benefits of working with large SDDs.

Previously, users had to choose between:

1. Working at low beam currents to maintain accuracy
2. Working with high beam current to get high count rates, risking beam damage and contamination

**Now X-Max offers users two more options:**


3. Work at low beam currents with high count rates and improved accuracy
4. Or decrease analysis times

**X-Max** The largest area SDD

# X-Max

## Size and productivity


Benefit immediately with increased counts from **X-Max**



These almandine garnet spectra were collected using a 10mm<sup>2</sup> SDD detector and the new 20mm<sup>2</sup>, 50mm<sup>2</sup>, and 80mm<sup>2</sup> X-Max detectors. The larger X-Max detectors collect more X-rays under the same collection conditions.

## True 100,000 cps performance

Many more analyses can be performed with 80mm<sup>2</sup> detectors at 100,000 cps



This logarithmic scale bar chart shows how high the beam current needs to be to achieve 100kcps on a traditional 10mm<sup>2</sup> SDD at 5kV and 20kV. It is a ten-fold increase over what is required to generate the same count rate on a new X-Max 80mm<sup>2</sup> SDD detector.

**SIZE**  
IT MATTERS  
IT MATTERS


# X-Max

## Large area SDD

### The effect of beam current on analytical data

#### Data collected with standard 10mm<sup>2</sup> SDDs


##### Insufficient data for usable nanoanalysis



##### Some detail but very few counts

Data collected with a 10mm<sup>2</sup> detector for 45 minutes at 3.3nA. Insufficient data has been collected which means a poor map that is not usable for analysis.

##### High beam current reduces accuracy




##### Lots of counts but no detail

Data collected with a 10mm<sup>2</sup> detector for 45 minutes at 36nA. The map shows lots of counts but the detailed information is blurred by a large spot size at this current.

##### Lots of counts and high quality data at low beam current

#### Data collected with 80mm<sup>2</sup> X-Max SDD



##### Lots of detail AND lots of counts



Data collected with an 80mm<sup>2</sup> detector for 45 minutes at 3.3nA. The mix map shows good spatial resolution and contains lots of data for good quantitative analysis.

# X-Max

## Size and accuracy

### Increasing accuracy by lowering accelerating voltage

X-Max offers unparalleled benefits when working at low kV (and low count rates).




### Maximise counts for accurate nanoanalysis with X-Max


These two maps were collected from the same area under the same conditions, changing only accelerating voltage. While superficially similar, points A, B, and C are actually very different. This is because reducing accelerating voltage (kV) reduces the interaction volume from which the X-rays are generated. As a result, the 20kV map shows X-ray signals from structures below the surface not seen in the 5kV map. Analysing at low kV benefits nanostructure analysis as it means that X-rays are only collected from the surface structures.

### Larger detectors cause less surface damage

#### Beam damage on sensitive samples




This image is of a polymer sample, taken at 0.2nA. Features on the surface and striation from sample preparation can clearly be seen.



The same sample imaged at 2nA shows extensive surface damage and distortion as a result of the higher beam current.

### Beam contamination



This image is of a sample from a nickel spinner bowl. Two areas have been examined, the one on the left using a large area detector, the one on the right with a traditional 10mm<sup>2</sup> detector. To collect the same data the sample had to be exposed for ten times longer with the smaller detector and as a result has far worse contamination.


# X-Max

## Technology leadership




### Unique large area sensors

- Discrete PentaFET® design offers the following benefits:
  - Protects FET from X-ray damage\*
  - Superior thermal design and FET gain stabilisation means detectors cool quickly and stabilise fast
  - Radial design for minimal ballistic deficit\*
- Polycrystalline Si construction process for extreme low leakage



### Unique magnetic trap design

- Allows large area sensor very close to samples
- Minimum sample to sensor distance = maximum capture angle
  - Typically 10 times greater than standard 10mm<sup>2</sup> SDDs
- Uses optimised rare earth element design




# X-Max

Pushing the boundaries



## Vacuum enclosed sensor

- No gas fill to absorb light elements
  - Nitrogen gas in other SDDs absorbs oxygen and other low energy X-rays
  - Vacuum enclosure increases analytical accuracy
- Long lifetime design
  - Uses OI vacuum technology
- Vacuum allows for greater cooling
  - Dramatically reduces thermal noise in large area SDDs
  - Not possible with gas-filled systems






The oxygen X-rays are absorbed by the nitrogen molecules present in the 'vacuum' so only a fraction of the X-ray signal gets through to the sensor.

The oxygen X-rays are not absorbed as there is no nitrogen present so you can see the true shape and height of the oxygen peak.

# X-Max

## Applications

### Low energy mapping of nano-particles





**X-Max** achieves excellent resolution. It also allows the user to do extraordinary nanoanalysis under conditions where other detectors struggle to produce sensible results. This means **X-Max** not only allows detection of small precipitates, but can map them as well.

Here we see a sample containing very small boride and carbide precipitates. The map clearly distinguishes between the different phases and even identifies the chemical composition of precipitates as small as 30nm. The spectrum above is from the highlighted 30nm carbide precipitate and includes data equivalent to a 0.5 second acquisition.

### Very high count rate maps

Texture analysis of fracture surface of a granite sample. This map, taken at a count rate of 400,000cps with an **X-Max** 80mm<sup>2</sup> detector, shows clearly how it is possible to map the different chemical compositions of a rough surface to see how the textured surface and chemical composition of the surface relate to each other.



# X-Max


## Superb performance

### Why does an **X-Max** EDS system give such exemplary performance?

- A great detector
  - Maximum real active area
  - Excellent resolution at high and low energies
  - Good throughput
  - ISO15632:2002 resolution compliance
- Excellent electronics
  - Highly linear analogue to digital conversion
  - Good rejection of simultaneous events (Pulse-pile-up)
  - Excellent signal to noise ratio
  - < 1eV change in peak resolution and position between 1 and 100,000cps
- Unique OI-owned interpretation algorithms
  - Complete - and correct - database of X-ray lines and edges
  - Software correction of sum-peaks
  - Highly accurate deconvolution, AutoID, and Quantitation

### What the system doesn't have:

- Shelf, tail, and shift correction - only real data is analysed
- Fan cooling - no vibration
- Forbidden elements - all elements from Be to Pu can be analysed
- Identification of absent elements - no false positives



At 80mm<sup>2</sup>, **X-Max** gives MnK $\alpha$  resolution that other SDDs can't match. Here we see an **X-Max** detector that achieves 123eV at MnK $\alpha$ ; better than the guaranteed specification of 124eV (Premium detector).

(Spectrum courtesy of P Russell, Appalachian State University. Recorded at 20kV, 168pA).

## X-Max

For **accuracy** / For **nano-analysis** / For **throughput** / For **productivity**

**SIZE**  
**IT MATTERS**  
ALL MATTERS

# X-Max

## In practice



Detector area	20mm <sup>2</sup>	50mm <sup>2</sup>	80mm <sup>2</sup>
<b>Resolution (eV):</b>			
MnK $\alpha$ resolution (Standard)	127	127	127
MnK $\alpha$ resolution (Premium)	124	124	124
<b>Operating angle</b>	0° to 45°		
<b>Temperature range</b>	10°C to 30°C		
<b>Altitude</b>	Sea level to 1,500m		

## Comparison of detector type

	Si(Li)	Other SDD	X-Max
<b>LN<sub>2</sub>-free</b>	NO	YES	YES
<b>Maximum single sensor active area</b>	30mm <sup>2</sup>	30mm <sup>2</sup>	80mm <sup>2</sup>
<i>Typical analysis conditions for most productive count rate</i>			
<b>Spot size</b>	Small	Large	Small
<b>Beam current</b>	Low	High	Low
<b>Accuracy at low kV</b>	High	Low	High
<b>Maximum throughput</b>	Low	High	High
<b>Productivity at low beam currents</b>	Low	Low	High

visit [www.oxford-instruments.com](http://www.oxford-instruments.com) for more information

This publication is the copyright of Oxford Instruments plc and provides outline information only, which (unless agreed by the company in writing) may not be used, applied or reproduced for any purpose or form part of any order or contract or regarded as the representation relating to the products or services concerned. Oxford Instruments' policy is one of continued improvement. The company reserves the right to alter, without notice the specification, design or conditions of supply of any product or service. Oxford Instruments acknowledges all trademarks and registrations. © Oxford Instruments plc, 2011. All rights reserved. Part no: OINA/160/C/1011



*The Business of Science®*

