

## X-MET<sup>®</sup> for heavy metals

Fast, reliable compliance screening!


### Prove compliance to industry regulations!

X-ray fluorescence (XRF) analyzers deliver fast, non-destructive screening for the detection of lead and other heavy metals restricted in toys, jewellery and consumer products. In seconds, prove compliance and gain peace of mind that only safe product will reach your customers.



#### X-MET for fast analysis of:

- Toys and other plastic goods for lead and other heavy metals (Pb, Cd, As, Sb, Ba, Cr, Se, Ni)
- PVC plastic: indication of phthalates
- Jewellery and metallic goods: Pb and Ni
- Paints and plastic raw materials
- RoHS compliance: Pb, Cd, Hg, total Cr and Br
- Packaging material: Cd, Pb and Hg



Silicon Drift Detector technology improves productivity!


Top of the range **X-MET5100** combines Oxford Instruments' groundbreaking Silicon Drift Detector (SDD) with a powerful 45kV X-ray tube. This cutting edge technology delivers a five times faster measurement speed, much better detection limits and significant accuracy improvement over conventional systems. Isn't it time you used **X-MET** to improve your productivity and screening confidence?



The Business of Science<sup>®</sup>

## Powerful user programmable interface for reliable testing

- Fast information for reliable Go/No-Go decisions
- User programmable elements, alarm limits and testing criteria
- Option to show only elements of interest or all elements
- High speed averaging – calculate averages of 2 – 50 measurements and save both individual and average results
- Save an unlimited number of results and spectra



Elem	ppm	STD	Alarm
Pb	1000	12	FAIL
Cd	235	0	FAIL
Hg	10	15	PASS
Sb	10	24	PASS
Se	10	5	PASS
Ba	1312	266	FAIL
As	10	6	PASS

## Flexibility with choice of analysis modes

- Auto detect mode will identify material type and automatically choose the correct analysis method:

**Plastics:** Fundamental Parameter (FP)  
Calculations for all common heavy elements in plastics: Pb, Cd, Hg, Br, Cr, Cl, Sr, Ni, Sb, As, Se, Ba etc.

**Metals:** Fundamental Parameter (FP)  
Calculations for alloy identification and analysis of all common heavy elements in metal alloys: Pb, Cd, Ni etc.

- Customer specific calibrations and optional software package for method development



## Easy and reliable

- Short learning curve
- User interface in >10 languages
- Easy data storage and reporting
- PDA based technology for flexibility
- CE, cCSAus certified

## Bench-Top Stand

- Improved precision
- Lower detection limits
- Optimized accuracy
- More comfortable operation



## Oxford Instruments Industrial Analysis

For more information please email:  
[industrial@oxinst.com](mailto:industrial@oxinst.com)

### UK

High Wycombe  
Tel: +44 (0) 1494 442255

### China

Shanghai  
Tel: +86 21 6132 9688

### Finland

Espoo  
Tel: +358 9 329 411

### Germany

Uedem  
Tel: +49 (0) 2825 93 83 -0

### Latin America

Tarpon Spring FL  
Tel: +1 978 369 9933 Ext. 220

### Singapore

Tel: +65 6337 6848

### North America

Concord MA  
TOLLFREE: +1 800 447 4717  
Tel: +1 978 369 9933

[www.oxford-instruments.com](http://www.oxford-instruments.com)

visit [www.oxford-instruments.com](http://www.oxford-instruments.com) for more information

This publication is the copyright of Oxford Instruments plc and provides outline information only, which (unless agreed by the company in writing) may not be used, applied or reproduced for any purpose or form part of any order or contract or regarded as the representation relating to the products or services concerned. Oxford Instruments' policy is one of continued improvement. The company reserves the right to alter, without notice the specification, design or conditions of supply of any product or service. Oxford Instruments acknowledges all trademarks and registrations. © Oxford Instruments plc, 2009. All rights reserved. Part no: OIIA/068/B/1009



*The Business of Science®*