

The search for an OES metals analyzer starts with the question...

...which instrument is best suited for my requirements?

...the search for the right supplier starts with Oxford Instruments!

Oxford Instruments, global leader in OES (Optical Emission Spectrometry) Arc/Spark technology, offers the widest range of OES analyzers.

With the latest in CCD technology, our instruments cover both ferrous and non-ferrous metal analysis, offering you flexibility in the element ranges required.

The two benchtop systems we recommend to the metals industry are:

FOUNDRY-MASTER UV

Fast, precise chemical analysis of all types of metal alloys. Ideally suited to foundries and metal producers.



FOUNDRY-MASTER COMPACT

Specifically designed for non-ferrous applications is ideal for routine analysis of metals and alloys.



Oxford Instruments Industrial Analysis

UK

Halifax Road, High Wycombe
Bucks, HP12 3SE England
Tel: +44 (0) 1494 442255
Fax: +44 (0) 1494 461033
Email: analytical@oxinst.com

China

Beijing
Tel: (8610) 6518 8160/1/2
Fax: (8610) 6518 8155
Email: info@oichina.cn

Finland

Espoo
Tel: +358 9 329 411
Fax: +358 9 3294 1300
Email: FI-Espoo_Info@oxinst.com

France

Saclay, Cedex
Tel: +33 (0) 1 69 85 25 24
Fax: +33 (0) 1 69 41 86 80
Email: analytical-info@oxford-instruments.fr

Germany

Uedem
Tel: +49 (0) 2825 93 83 -0
Fax: +49 (0) 2825 93 83 -100
Email: DE-Uedem_Info@oxinst.com

Japan

Tokyo
Tel: +81 (0) 3 5245 3591
Fax: +81 (0) 3 5245 4466/4477
Email: oikkma@oxinst.co.jp

Latin America

Clearwater FL
Tel: +1 727 538 7702
Fax: +1 727 538 4205
Email: oxford@gate.net

Singapore

Tel: +65 6337 6848
Fax: +65 6337 6286
Email: asiasales@oxinst.com

North America

Elk Grove Village IL
Tel: +1 847 439 4404
Fax: +1 847 439 4425
Email: sales@rmys.oxinst.com

www.oxford-instruments.com

Click onto www.oxford-instruments.com for more information

Oxford Instruments, at High Wycombe, UK, operates Quality Management Systems approved to the requirements of BS EN ISO 9001. This publication is the copyright of Oxford Instruments Analytical Limited and provides online information only which (unless agreed by the company in writing) may not be used, copied or reproduced for any purpose or form part of any order or contract or be regarded as a representation relating to the products or services concerned. Oxford Instruments' policy is one of continued improvement. The company reserves the right to alter, without notice, the specification, design or conditions of supply of any product or service. Oxford Instruments acknowledges all trade marks and registrations.

© Oxford Instruments Analytical Ltd, 2008. All rights reserved.



Certificate No FM29142

